

## SHAKER-INCUBATORS For 2 or 4 Microplate



### Model: MOB-02A/04A, 2/4 Thermo Shaker Incubator for plates

MOB-02A/MOB-04A are designed for 1–4 standard 96 well plates in the thermostat mode. MOB-02A/MOB-04A have incubation and shaking functions.

#### Features:

- Heating plates above and below the microplates
- Exact correspondence of the designed & actual temp.
- LCD display. It is easy to set up and use.

Model	MOB-02A	MOB-04A
Temperature range	RT+5°C~70°C	
Timer range	1min-99h59min	
Temp. accuracy	±0.5°C	
Display accuracy	±0.1°C	
Shaking speed	100-1200 rpm	100-1500 rpm
Mixing Orbit	2mm	
Heating Time (20 to 70°C)	≤20min	
Platform capacity	2 microtitre plates	4 microtitre plates
Dimensions	280x270x150mm	320x350x185mm
Weight	6.5kg	9kg