



**AMP USR *compact*** - Urine Strip Reader



Reliable

# Urine Diagnostics

Simple and Efficient



## AMP Urine Strips

Performing screening tests in hospitals, health care centres and general practices with AMP Urine Strips are a simple, inexpensive and efficient tool to identify risk patients requiring further diagnostic examination. The test strips can be evaluated both visually and by using adequate automated instruments.

Use of only the best raw materials and an efficient quality control system ensure AMP Urine Strips to report results of optimal accuracy for each parameter. The arrangement of the individual test pads ensures the efficient prevention of cross contamination.

### AMP Urine Test Strips and Quality Control Material

|        | URO | BIL | KET | BLD | PRO | NIT | LEU | GLU | SG | pH | ASA | MAL | B | CRE | CA |
|--------|-----|-----|-----|-----|-----|-----|-----|-----|----|----|-----|-----|---|-----|----|
| US-1   |     |     |     |     |     |     |     | •   |    |    |     |     |   |     |    |
| US-1P  |     |     |     |     | •   |     |     |     |    |    |     |     |   |     |    |
| US-1K  |     |     | •   |     |     |     |     |     |    |    |     |     |   |     |    |
| US-1B  |     |     |     | •   |     |     |     |     |    |    |     |     |   |     |    |
| US-2   |     |     |     |     | •   |     |     | •   |    |    |     |     |   |     |    |
| US-2K  |     |     | •   |     |     |     |     | •   |    |    |     |     |   |     |    |
| US-2N  |     |     |     |     |     | •   | •   |     |    |    |     |     |   |     |    |
| US-2M  |     |     |     |     |     |     |     | •   |    |    |     |     | • |     |    |
| US-3   |     |     |     |     | •   |     |     | •   |    | •  |     |     |   |     |    |
| US-3K  |     |     | •   |     | •   |     |     | •   |    |    |     |     |   |     |    |
| US-4   |     |     |     |     | •   |     |     | •   | •  | •  |     |     |   |     |    |
| US-8   | •   | •   | •   | •   | •   | •   | •   | •   |    | •  |     |     |   |     |    |
| US-8KS |     |     | •   | •   | •   | •   | •   | •   | •  | •  |     |     |   |     |    |
| US-8US | •   |     |     | •   | •   | •   | •   | •   | •  | •  |     |     |   |     |    |
| US-8KU | •   | •   | •   |     | •   | •   | •   | •   |    | •  |     |     |   |     |    |
| US-10  | •   | •   | •   | •   | •   | •   | •   | •   | •  | •  |     |     |   |     |    |
| US-11  | •   | •   | •   | •   | •   | •   | •   | •   | •  | •  | •   |     |   |     |    |
| US-11M | •   | •   | •   | •   | •   | •   | •   | •   | •  | •  | •   |     | • |     |    |
| US-12  | •   | •   | •   | •   | •   | •   | •   | •   | •  | •  | •   | •   |   |     |    |
| US-13  | •   | •   | •   | •   | •   | •   | •   | •   | •  | •  | •   | •   | • | •   |    |
| US-14  | •   | •   | •   | •   | •   | •   | •   | •   | •  | •  | •   | •   | • | •   | •  |

**AMP UroTrol:**

Quality Control Material available in 2 levels (negative and positive)

### AMEDA Labordiagnostik GmbH

Krenngasse 12

8010 Graz, Austria

Phone +43 316 69 80 69

Fax +43 316 69 80 69 12

E-mail office.graz@amp-med.com

[www.amp-diagnostics.com](http://www.amp-diagnostics.com)



AMP Polska sp. z o.o.  
AMP Diagnostics s.r.l.

Bolesława Prusa 37/1, 30-117 Kraków, Poland, Phone: +48 12 294 02 78, E-mail: office.krakow@amp-med.com  
Str. Fabricii 118-A, 400632 Cluj-Napoca, Romania, Phone: +40 264 425 934, E-mail: office.cluj@amp-med.com